

Not **All** Feeding Tubes Are Created Equal!



CORFLO[®]

CORFLO[®] Enteral Feeding Tubes from VIASYS MedSystems offer advantages to your patients **unlike any other feeding tubes available**, including:

● Patented **ULTRA** anti-clog outlet port significantly reduces the occurrence of clogging.

● The largest selection of adult, pediatric, and neonatal feeding tubes to allow one-stop sourcing for all your feeding tubes hospital or system-wide.

● Patented **TRAK** feeding tubes: The **ONLY** feeding tubes that can be used with the **CORTRAK[™]** Enteral Access System.



The CORFLO® Advantage

The Ultimate in Patient Safety:

- All CORFLO feeding tubes are printed with centimeter marks for increased patient safety and intubation accuracy.
- The entire CORFLO polyurethane feeding tube line is color coded to correspond to varying lengths.
- Stylet lengths are matched to the CORFLO feeding tubes.

Clinician Convenience and Patient Comfort:

- Flow-through stylets to flush, auscultate and aspirate during the insertion procedure
- Water activated lubricant on all polyurethane feeding tubes for greater patient comfort during intubation and effortless stylet removal
- 100% radiopaque for X-ray confirmation
- Weighted tubes use multiple tungsten weights allowing for flexibility and greater comfort during intubation
- "Y" ports on all CORFLO polyurethane feeding tubes minimize touch contamination and makes flushing easy

VIASYS MedSystems feeding tubes are the standard by which all other feeding tubes are measured!

VIASYS MEDSYSTEMS 20-9431 10 Fr 43"-109cm 100 + + + + 95 +

FR Size/ Length	5 FR STERILE	6 FR	8 FR	10 FR	12 FR
15" (38cm)	20-1155* 20-PVC155*	20-1156*	20-PVC158		
22" (56cm)	20-9225* 20-C1225*	20-9226 20-1226 20-8226 20-C1226	20-9228 20-C1228		
36" (91cm)	20-1365*	20-9366 20-1366 20-8366 20-7366	20-9368 20-1368 20-8368 20-7368	20-9361 20-1361 20-7361 20-5361	20-1362 20-5362
43" (109cm)			20-9438TRAK 20-9438 20-1438 20-7438 20-5438 20-6438 20-4438 20-3438 20-2438 20-0438	20-9431TRAK 20-9431 20-1431 20-7431 20-5431 20-6431 20-4431 20-3431 20-0431	20-9432 20-1432 20-7432 20-5432 20-4432 20-0432
55" (140cm)	All Feeding Tubes packaged 10/case, except Infant PVC (50/case)		20-9558	20-9551 20-9551TRAK 20-3551 20-2551	

CORFLO®-ULTRA NG:** designed for intubations where anti-clogging and ease of insertion are critical.
CORFLO-CONTROLLER NG: developed for maximum control during the intubation procedure especially with Endoscopy and Fluoroscopy placements.
CORFLO-INFANT PVC NG: manufactured from premium polyvinyl chloride, packaged sterile for short-term feeding for infants.

CORFLO-ULTRA Lite NG Non-wt. w/stylet
CORFLO-ULTRA Lite NG Non-wt. w/o stylet
CORFLO-ULTRA Lite Clear NG STERILE, Clear, non-wt. w/o stylet

CORFLO-ULTRA Pedi NG 1 g & 1.5g Pediatric wt. w/stylet
CORFLO-ULTRA NG 3 g wt. w/stylet
CORFLO-ULTRA NG 3 g wt. w/o stylet
CORFLO-ULTRA PIII NG 3 g wt. w/stylet
CORFLO-ULTRA 7 NG 7 g wt. w/stylet
CORFLO-CONTROLLER NG 3 g wt. w/stylet
CORFLO-CONTROLLER PIII NG 3 g wt. w/stylet
CORFLO-CONTROLLER 7 NG 7 g wt. w/stylet
CORFLO-INFANT P NG STERILE, Plastic, non-wt. w/o stylet

CORTRAK ENTERAL ACCESS SYSTEM NG Non-wt. w/stylet

*STERILE
 **All CORFLO-Ultra feeding tubes have our unique ANTI-CLOG feeding port.
 Boluses: approximate weights

Enterally speaking, nobody delivers better!®

VIASYS MedSystems
 100 Chaddick Dr.
 Wheeling, IL 60090

1-800-323-6305
 1-847-403-3400 • FAX 1-847-541-9526
www.viasyshealthcare.com

VIASYS
 HEALTHCARE
 MEDSYSTEMS

CE
 ISO 13485
 5841D 09/15