

SURE SEAL® face masks incorporate specially designed nose area and comfortable, gas-tight, low pressure seal

air cushion masks



SURE SEAL® Masks

- Colour-coded retaining rings are removable for hand-held procedures
- Available with optional inflation port for customised fit and comfort
- Individually packaged 20 to a case

Cat No	Size	Connector Size (mm)	Inflation Port	Ring Colour
HUD1275	Medium Adult	22		Red
HUD1276	Large Adult	22		Blue
HUD1278	Infant	15	•	Pink
HUD1281	Medium Adult	22	•	Red
HUD1282	Large Adult	22	•	Blue

• masks cannulas + tubing • nebulisers • HME/filters • **airway management** • adaptors + miscellaneous

HUDSON RCI®



LIFESAVER® reusable resuscitators

The LIFESAVER® Reusable Resuscitator is designed to meet ISO and ASTM standards

● masks cannulas + tubing ● nebulisers ● HME/filters ● **airway management** ● adaptors + miscellaneous



The LIFESAVER® Reusable Resuscitator is designed to meet ISO and ASTM standards.

- Autoclavable (except reservoir bag)
- Easy to assemble
- Collapsible, for compact, convenient storage
- Non-rebreathing valve offers low resistance to both inspiratory and expiratory gas flows
- Cat No HUD1669 adaptor positioned between mask and bag, may be used for pressure or CO₂ monitoring
- Sold individually

adult resuscitator

Cat No HUD5340

ADULT BAG

- Mask and reservoir bag not included

Cat No HUD5345

ADULT BAG

- Includes mask, reservoir bag, and reservoir valve assembly

child resuscitator

Cat No HUD5341

CHILD BAG

- Mask and reservoir bag not included

Cat No HUD5346

CHILD BAG

- Includes mask, reservoir bag, and reservoir valve assembly

replacement parts

Cat No	Description
HUD51007	Oxygen Reservoir, 2500ml

replacement masks

Cat No	Size
HUD5190	Small Infant
HUD5191	Large Infant
HUD5192	Small Child



HUDSON RCI®

The LIFESAVER® Disposable Manual Resuscitators combine quality, features and performance

LIFESAVER® disposable resuscitators



LIFESAVER® Disposable Manual Resuscitators combine quality, features and performance.

- Designed to meet ISO and ASTM standards
- Unique right-angle O₂ tubing connector swivels 360 degrees
- 100% latex-free
- Individually packaged 6 to a case

adult resuscitation bags

Available with a variety of options including a capped pressure monitoring port/60cm pop-off valve style.

- All configurations include swivel flow diverter for connection of PEEP valve
- Available with optional PEEP valve

Cat No	Mask (#1276)	Valved Mask (#1282)	PEEP Valve	Pop-Off Valve
HUD5372	•			
HUD5377	•			•

paediatric and infant resuscitation bags

All paediatric/child and infant/neonatal bags include a capped pressure monitoring port and 40cm pressure relief pop-off valve.

- All configurations include swivel flow diverter for optional connection of PEEP valve
- Available with optional PEEP valve

Cat No	Size	Paediatric Mask (#1273)	Neonate Mask (#1271)	PEEP Valve
HUD5362	Infant		•	
HUD5367	Paediatric	•		

disposable PEEP valve (30mm) – adult

Cat No HUD5385

- Adjustable 5 to 20cm H₂O
- ISO standard 30mm ID connector
- Attaches to any resuscitation bag flow diverter with 30mm connection
- Packaged 10 to a case

• masks cannulas + tubing • nebulisers • HME/filters • airway management • adaptors + miscellaneous

HUDSON RCI®



oral airways



colour-coded airways

- Semirigid, non-toxic polyethylene with colour-coded bite block
- Guedel-style enclosed centre channel
- Individually packaged 48 to a case

Cat No	No	Size	Colour
HUDI161	11	110mm	Orange
HUDI160	10	100mm	Red
HUDI159	9	90mm	Yellow
HUDI158	8	80mm	Green
HUDI157	7	70mm	White
HUDI156	6	60mm	Black
HUDI155	5	50mm	Blue
HUDI154	4	40mm	Pink



CATH-GUIDE® airways

- Flexible vinyl with a rigid bite block
- Guedel-style with three internal channels
- Individually packaged 48 to a case

Cat No	No	Size
HUDI165	12	120mm
HUDI166	11	110mm
HUDI167	10	100mm
HUDI164	9	90mm
HUDI168	8	80mm
HUDI169	7	70mm
HUDI170	6	60mm
HUDI171	5.5	55mm



PocketPak™ resuscitator

Cat No HUDI 187

The PocketPak™ emergency mask kit provides effective protection for mouth-to-mouth resuscitation and fits standard manual resuscitators. The collapsible mask is easy to use and convenient to carry.

PocketPak™ features:

- Mask with pliable, clear dome to provide secure seal
- Capped port for delivery of supplemental oxygen
- Low-resistance valve and particle filter
- Foldable extension tube
- Packaged in a sturdy 'clam shell' case,
- 4 masks to a convenient storage box



BiteGard™

Cat No HUDI 140

Disposable, compressible molar bite block provides a cost-effective solution for patients who clench their teeth.

- Prevents damage to laryngeal mask airways (LMA), bronchoscopes and other valuable medical equipment
- Individually packaged 50 to a case



HUDSON RCI®

