

# oxygen masks



## medium concentration masks

### Cat No HUD1040

ADULT

- Clear, soft vinyl for patient comfort and visual assessment
- Adjustable nose clip assures comfortable fit
- Complete with 2.1m oxygen supply tubing

### Cat No HUD1035

PAEDIATRIC

- Same as Cat No HUD1040 except paediatric size mask

### Cat No HUD1048

ADULT

- Same as Cat No HUD1040 but without tubing

*Masks are individually packaged 50 to a case*



## medium concentration masks – elongated

### Cat No HUD1041

ADULT

- Under-the-chin design for excellent fit, wide range of face sizes
- Clear, soft vinyl for patient comfort and visual assessment
- Adjustable nose clip assures comfortable fit
- Complete with 2.1m oxygen supply tubing

### Cat No HUD1049

ADULT

- Same as Cat No HUD1041 but without tubing

### Cat No HUD1042

PAEDIATRIC

- Same as Cat No HUD1041 except paediatric size mask

### Cat No HUD1050

PAEDIATRIC

- Same as Cat No HUD1042 but without tubing

*Masks are individually packaged 50 to a case*

# oxygen masks



## high concentration masks – elongated

### Cat No HUD1007

ADULT

- Under-the-chin design for excellent fit on wide range of face sizes
- Clear, soft vinyl for patient comfort and visual assessment
- Adaptor swivels to accommodate patient position
- Adjustable nose clip assures comfortable fit
- 750ml reservoir bag
- Complete with 2.1m oxygen supply tubing

### Cat No HUD1011

PAEDIATRIC

- Same as Cat No HUD1007 except paediatric size mask

*Masks are individually packaged 50 to a case*



## non-rebreathing masks – with safety vent

### Cat No HUD1059

ADULT

- Low-resistance, latex-free, check valve prevents rebreathing and allows exhaled gas to escape
- Safety vent allows for entrainment of room air
- Adjustable nose clip assures comfortable fit
- 750ml reservoir bag
- Complete with 2.1m oxygen supply tubing

### Cat No HUD1069

ADULT

- Same as Cat No HUD1059 but without tubing

### Cat No HUD1058

PAEDIATRIC

- Same as Cat No HUD1059 except paediatric size mask

*Masks are individually packaged 50 to a case*

**HUDSON RCI®**



# oxygen masks

● masks cannulas + tubing ● nebulisers ● HME/filters ● airway management ● adaptors + miscellaneous



## non-rebreathing masks

**Cat No HUD1060**

ADULT

- Low-resistance, latex-free check valve prevents rebreathing and allows exhaled gas to escape
- Adjustable nose clip assures comfortable fit
- 750ml reservoir bag
- Complete with 2.1m oxygen supply tubing



## three-in-one masks

**Cat No HUD1061**

ADULT

- Three-in-one masks can be used as a medium concentration, high concentration, or non-rebreathing mask
- Includes mask with latex-free flapper valve, non-rebreathing bag assembly, 2.1m oxygen supply tubing, and oxygen connector
- Adjustable nose clip assures comfortable fit
- 750ml reservoir bag

*Masks are individually packaged 50 to a case*

**Cat No HUD1937**

ADULT

- Same as Cat No HUD1061 but with STAR LUMEN® tubing



**HUDSON RCI®**

# oxygen masks + air entrainment masks



**MULTI-VENT® Masks**

**Cat No HUD1088**

ADULT

- Safe, simple delivery of variable oxygen concentrations
- Each mask includes: colour-coded, low and medium concentration diluters
- Locking ring secures flow setting
- Includes adaptor for high humidity entrainment
- Complete with 2.1 m oxygen supply tubing

**Cat No HUD1089**

PAEDIATRIC

- Same as Cat No HUD1088 except paediatric size mask

*Masks are individually packaged 50 to a case*

*Dilutors are not sold separately*



**SELECT-A-VENT™ Masks**

- Safe, simple delivery of variable oxygen concentrations
- Features six colour-coded diluters: 24%, 28%, 31%, 35%, 40% and 50%
- Diluters not sold separately
- Includes adaptor for high humidity entrainment
- Complete with 2.1 m oxygen supply tubing

**Cat No HUD1098**

ADULT

- Kit includes mask, high humidity adaptor; oxygen tubing and six colour-coded diluters

**Cat No HUD1096**

PAEDIATRIC

- Same as Cat No HUD1098 except paediatric size mask

*Masks are individually packaged 50 to a case*

*Dilutors are not sold separately*

**HUDSON RCI®**



● masks cannulas + tubing ● nebulisers ● HME/filters ● airway management ● adaptors + miscellaneous

# aerosol masks



## aerosol masks – elongated

### Cat No HUD1083

#### ADULT

- Clear, soft vinyl for patient comfort
- Adjustable nose clip assures comfortable fit
- Specifically designed for aerosol therapy
- Does not include tubing
- Accepts 22mm ID aerosol tubing

### Cat No HUD1085

#### PAEDIATRIC

- Same as Cat No HUD1083 except paediatric size mask

*Masks are individually packaged 50 to a case*



## face tents

### Cat No HUD1095

- For high-humidity aerosol therapy
- Distinctive under-the-chin design made of clear, soft vinyl offers patient comfort and ease of speech
- Does not include tubing
- Accepts 22mm ID aerosol tubing

*Masks are individually packaged 50 to a case*

# aerosol masks + oxygen catheters + cannulas



**tracheostomy masks**

**Cat No HUD1075**

ADULT


- For tracheostomy and laryngectomy aerosol therapy
- Tubing connector swivels 360 degrees for ease of positioning
- 22mm OD connector accepts 22mm corrugated tubing and nebuliser tees

**Cat No HUD1076**

PAEDIATRIC

- Same as Cat No HUD1075 except paediatric size mask

*Masks are individually packaged 50 to a case*



**oxygen catheters**

**Cat No HUD1147**

14 FRENCH

- 40cm flexible oxygen catheter made of smooth, odourless material
- Apertures are spaced to provide optimal patient comfort
- Individually packaged 50 to a case

**Cat No HUD1145**

10 FRENCH

- Same as Cat No HUD1041 but 10 French

**HUDSON RCI®**



• masks cannulas + tubing • nebulisers • HME/filters • airway management • adaptors + miscellaneous

Hudson RCI offers an extensive line of nasal oxygen cannulas to accommodate a variety of patient needs in all care settings

# nasal cannulas

● masks cannulas + tubing ● nebulisers ● HME/filters ● airway management ● adaptors + miscellaneous



## over-the-ear cannulas

The popular Hudson RCI over-the-ear style cannula which has set the industry standard for decades is cost effective and offered in a wide range of adult configurations.

- Flexible, contoured lip tab provides a high level of stability and patient comfort

- Over-the-ear design provides comfortable and secure fit
- Offered with standard (non-flared) or flared nasal prongs
- Available with standard or universal connector and standard or crush-proof STAR LUMEN® tubing in a variety of lengths
- Cannulas individually packaged, case quantities are listed below

## over-the-ear cannulas configurations

Cat No	Standard Nasal Tip	Flared Nasal Tip	Tubing Length	STAR LUMEN® Tubing	Standard Tubing	Standard Connector	Universal Connector	Case Qty
HUD1109	•		none					50
HUD1103	•		2.1m		•	•		50
HUD1110	•		2.1m	•		•		50
HUD1108		•	none					50
HUD1104		•	2.1m		•	•		50
HUD1111		•	2.1m	•		•		50



**HUDSON RCI®**

Hudson RCI offers an extensive line of nasal oxygen cannulas to accommodate a variety of patient needs in all care settings

# nasal cannulas



## SOFTECH® nasal cannulas

The SOFTECH® nasal cannula is produced by a special manufacturing process to provide an exceptionally soft cannula.

- Extremely lightweight and flexible
- Anatomically designed to enhance patient comfort and acceptance
- Over-the-ear design optimises fit and stability while maximising patient comfort

- Available with standard or universal tubing connector which allows direct attachment to flow meter
- Most configurations include crush-proof and kink-proof tubing for continuous flow
- Available without tubing
- Cannulas individually packaged in case quantities listed below

## nasal cannulas configurations

Cat No	Size	Tubing Length	STAR LUMEN® Tubing	Standard Tubing	Standard Connector	Universal Connector	Case Qty
HUD1820	Adult	2.1m	•		•		50
HUD1821	Adult	none					50
HUD1826	Paediatric	2.1m	•			•	50
HUD1828	Infant	2.1m	•		•		50
HUD1838	Paediatric	2.1m	•		•		50

**HUDSON RCI®**



# infant CPAP cannulas

The Hudson RCI Infant Nasal Prong CPAP cannula system is designed to reduce trauma associated with the delivery of infant nasal CPAP

● masks cannulas + tubing ● nebulisers ● HME/filters ● airway management ● adaptors + miscellaneous



### Infant Nasal Prong CPAP Cannula

- Soft, anatomically curved prongs enhance fit and minimize nasal septal necrosis
- Luer fitting on expiratory connector allows proximal airway pressure monitoring
- Six prong sizes allow greater choice for appropriate sizing for each infant
- Available in a 'set' for use with most mechanical ventilators
- Educational videotape 'Infant Nasal Prong CPAP' available
- Individually packaged 10 to a case

### CPAP sets

Sets are intended for use with mechanical ventilators and include:

- Cannula
- Knit cap
- Inspiratory elbow connector
- 2 x 15.2cm Velcro® sections
- Expiratory elbow connector
- 2 x 10mm to 7.5mm adaptors

Infant Weight Range	Suggested Cannula Size	French Size	Set Cat No
<700 grams	0	10	HUD1690
700 - 1,250 grams	1	12	HUD1691
1,250 - 2,000 grams	2	13	HUD1692
2,000 - 3,000 grams	3	15	HUD1693
>3,000 grams	4	16	HUD1694
1 - 2 years	5	20	HUD1695



**HUDSON RCI®**

# oxygen supply tubing

• masks cannulas + tubing • nebulisers • HME/filters • airway management • adaptors + miscellaneous



**oxygen tubing**

Hudson RCI oxygen supply tubing is available in a variety of configurations to meet the needs of patients and clinicians in all care settings.

- Available with standard or universal connector which attaches directly on to flowmeter eliminating need for an adaptor
- Distinctively designed STAR LUMEN® tubing features a channelled inner lumen to resist kink and occlusion

**oxygen tubing configurations**

Cat No	Standard Tubing	STAR LUMEN® Tubing	Standard Connector	Universal Connector	Tubing Length (m)	Case Qty
HUD1112		•	•		2.1	50
HUD1115	•		•		2.1	50
HUD1118		•	•		4.3	50
HUD1119		•	•		7.6	25

**HUDSON RCI®**

