

small volume nebulisers

MICRO MIST® nebuliser

The MICRO MIST® small volume nebuliser is designed for performance and economy and includes the design features needed by today's respiratory care practitioners. It can be used for hand-held or in-line treatments.

- Optimal particle size delivered at efficient nebulisation rates (MMAD = 2.1 microns)
- Consistent performance at angles up to 90 degrees allowing medication delivery to patients positioned horizontally or with special requirements
- Jet stays in place unless intentionally removed
- Easy-seal, threaded cap and 8cc capacity anti-spill jar to prevent leakage

- Available in NEB-U-MASK® configuration with choice of adult or paediatric mask
- Available with standard or universal tubing connector
- Individually packaged 50 to a case

Cat No HUD1888

MICRO MIST NEBULISER KIT, 2+2

Convenient kit, ideal for home care applications.

- Each kit contains:
 - 2 x Cat No 1880
 - 2 x Cat No 1883
- Packaged complete, 25 kits to a case



nebuliser configurations

Cat No	Nebuliser	Tee	Mouthpiece	Reservoir Tube	2.1m Tubing	Standard Connector	Universal Connector	Adult Mask (#1083)	Paed Mask (#1085)
HUD1880	•								
HUD1881	•	•	•						
HUD1882	•	•	•		•	•			
HUD1883	•	•	•	•	•	•			
HUD1885	•				•	•		•	
HUD1886	•				•	•			•



HUDSON RCI®

small volume nebulisers



ISO-NEB™ Filtered Nebuliser System

Cat No HUD1755

- Designed to reduce environmental contamination during aerosol therapy
- Inspiratory and expiratory flow are isolated by latex-free one-way valves with exhalation directed through an efficient filter system
- Electrostatically charged pathogen filter is proven* 99.999% bacterial / 99.99% viral filtration efficient
- Anti-spill UP-DRAFT II OPTI-NEB™ nebuliser produces fine, dense particles for specific site therapy
- High-output, 8cc nebuliser operates vertically or at a 45° angle
- Particle sizes range from 1 to 2 microns mean diameter* for optimum deposition
- Individually packaged 20 to a case

**Verified by independent laboratory tests*



CORR-A-FLEX® Tubing

Cat No HUD1680

- Continuous 100ft roll of 22mm disposable corrugated tubing with cuttable 6 inch sections
- Packaged in protective, convenient, plastic-lined dispenser box
- Flexible, lightweight tubing
- Tapered cuffs securely attach to most 22mm respiratory therapy applications

HUDSON RCI®



● masks cannulas + tubing ● **nebulisers** ● HME/filters ● airway management ● adaptors + miscellaneous

large volume nebulisers

masks cannulas + tubing • nebulisers • HME/filters • airway management • adaptors + miscellaneous



AQUAPAK® prefilled nebulisers

Versatile, sterile systems to provide cold or heated nebulisation. Choose from a wide variety of reservoirs, adaptors and solutions to meet specific patient requirements efficiently and economically.

- Prefilled sterile reservoirs come in four convenient sizes; 440ml, 760ml, 1070ml and 2200ml
- Reservoirs may be used with other Hudson RCI humidification and nebulisation products, making the system truly versatile

nebuliser adaptors

Nebuliser adaptors designed for use with prefilled sterile reservoirs

Cat No HUD031-28 B

028 ADAPTOR (YELLOW)

- Venturi-style adaptor allows the selection of oxygen concentrations from 28% to 98%

AQUAPAK® prefilled reservoir and nebuliser adaptor configurations

Cat No	440 ml	760 ml	1070 ml	2200 ml	Sterile Water	0.45% Saline	0.90% Saline	028 Adaptor	033 Adaptor	4 per case	10 per case	20 per case	50 per case
HUD004-00	•				•							•	
HUD037-00		•			•						•		
HUD031-28B								•					•



HUDSON RCI®