

Infu-Stat™ Pressure Infuser

- 100% Latex Free
- Cost Competitive
- User Friendly
- One hand operation



Infu-Stat™ Pressure Infuser



- Latex Free
- Loads easily
- One hand operation for inflation or deflation
- Large bore stopcock for rapid inflation and deflation
- Color-coded gauge is easy to read and maintains reliable pressure
- The safety valve eliminates the possibility of over inflation and releases air at 325mm/Hg
- Sizing is color-coded for easy identification
- Graphics are sharp and easy to read



Item #	Size	Units
MTY 305	500ml	5/box
MTY 310	1000ml	5/box
MTY 330	3000ml	2/box

